TIPS:

Your suggestions and comments for COSTWAY are really important to us!
We sincerely solicit you to go back to our shop and leave a good rating in just a simple click. It would be quite encouraging if you could kindly do so like below:

★★★★★ Great products so far

February 24, 2018

Great products so far. Fast delivery, easy setup, and working without any issues.

With your inspiring rating, COSTWAY will be more consistent to offer you EASY SHOPPING EXPERIENCE, GOOD PRODUCTS and EFFICIENT SERVICE!

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2 Section Portable Massage Table

Precautionary Instructions

1. Maintain the table periodically. Make sure all screws are secure and the cable has not deteriorated or frayed. Make sure there is no cracks on the aluminum structure.
2. Do not overload the table beyond the recommended weight limit. This could result in serious injury to the user. This products is tested to hold a max working load of 250 kgs.
3. Do not use alcohol based cleaners when cleaning the vinyl. This could cause damage to the vinyl material. Use warm soapy water and wipe dry.
4. Make sure when using the table that all four legs are at the same level. Using the legs at different levels can cause the table to be unstable and cause injury.
5. Do not allow children to assemble or disassemble the massage table for you unless under supervision.
6. Be careful when using sharp objects as they can pierce the covering. Keep the table's covering away from being contacted by the sharp items.

1. Handle 2. Copper button 3. Hinge 4. Hanged armrest
 5. Headrest 6. Cable 7. Side arm supports

Setting Up Your Table

Steps are as below:

1. Release the buckle.
2. Lay the table on its side with four rubber feet on the floor, open up the table with four legs extend as full as possible.
3. Take out all accessories and then lift the table up until all legs are on the floor.
4. Release the side knob of the headrest cradle and then press and fix it once the headrest is adjusted to the required position..
5. Put on the hanged armrest or sling on the headrest cradle, Put the face cushion on the cradle, and then insert the two side arm supports into the two holes beside the table.
6. Now the table is ready for your use.

Re: The height of table is adjustable by pressing the copper buttons on the leg. Four legs should be done in the same way one by one and make sure all legs are adjusted into same height level.

Close and transport table

1. Take out all accessories..
2. Lay the table on its side with four rubber feet on the ground.
3. Put all accessories under the elastic strap under the top veneer.
4. Fold the table up and lock the buckle.
5. Put the massage table into the carry bag.And Carry the table like suitcase by grasping the handle.

Warning: If the frame/legs/table structure is damaged in anyway, then you should not use this table. Using this table when damaged may cause harm or serious injury to you or your patient. You should seek repair immediately.

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DE office:FDS GmbH,Neuer Höltingbaum 36,22143 Hamburg,Deutschland